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The informations here contained are subject to modifications without warning.

No part of this document could be photocopied, reproduced or translated in other language without the written consent of LED SpA.

In case of necessity of technical assistance contact your own retailer

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WARNINGS	25
INTRODUCTION	27
GENERAL INFORMATION	30
CONTROLS	32
<i>TECHNICAL CHARACTERISTICS</i>	<i>33</i>

<i>SIZE</i>	33
<i>STANDARD COMPOSITION</i>	33
COUNTER-INDICATIONS AND COLLATERAL EFFECTS	34
OPERATIONAL INSTRUCTIONS	35
TRANSDUCERS PLACEMENT	37
MAINTENANCE	38
APPLICATIONS	39
THERAPEUTIC INDICATIVE PROTOCOLS	40

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WARNINGS

It is very important that this manual of instructions is preserved together the unit for any future consultation. If the unit must be sold or transferred to other consumer, make sure that the manual comes furnished together, so that the new consumer could be aware of the operation of the unit and of the relative instructions.

THESE INSTRUCTIONS HAS STAYED COMPILED FOR YOUR SAFETY AND FOR THAT OF THE OTHERS, THEREFORE WE ASK YOU TO READ THEM ATTENTIVELY BEFORE TO INSTALL AND TO USE THE UNIT.

- After have removed the packing make sure of the integrity of the unit, in case of doubt don't use the unit and address to personal professionally qualified.
- The electric safety of this apparatus is insured only when the same is correctly connected to an efficient net for the connection to the earth in conformity with the actual safety requirements. It is necessary to verify this fundamental safety requisite and, in case of doubt, to require an accurate control of the plant from part of qualified personnel. The manufacturer cannot be considered responsible for possible damages caused from the lack of efficient connection to earth of the installation.
- Before connect the unit ascertain that the label data (on the back panel of the unit) are in conformity with those of the net of electric distribution.
- In case of incompatibility between the mains socket and the power cable of the unit replace it with an other of proper type. In general it is not advised the use of adapters, multiple connectors and/or extension cables. Should their use becomes essential, it is necessary to use simple or multiples adapters conforming to the actual safety norms only.
- The use of an any electric apparatus involves the observance of any fundamental rules. Particularly:
 - don't touch the apparatus with hands excessively wet
 - doesn't let exposed the apparatus to atmospheric agents (rain, sun, etc.)

- Don't leave the apparatus uselessly inserted. Switch off the general interrupter when the apparatus has not used.
- This unit will be destined only to the use for which it has been expressly conceived. Each other use has to be considered improper and therefore dangerous. The manufacturer cannot be considered responsible for possible damages due to improper, unreasonable or wrong uses.
- It is dangerous to modify or try to modify the characteristics of this unit.
- Before effect any operation of cleaning or maintenance, disconnect the unit from the electric net, or disconnecting it from the wall or switching off the general interrupter of the plant.
- In case of breakdown or bad operation of the unit, switch off it. For the possible reparation address only to a center of authorized assistance and ask for the use of original spare parts. Missing the fulfillment of the previous rules could compromise the safety of the unit and of the consumer, beyond to make to decay each form of guarantee.
- The use of this apparatus is not recommended in patients:
 - carriers of pace-maker or of other electronic installed device, or in proximity of such carrier patient
 - with serious unbalances of the arterial pressure
 - with serious illness of the nervous system
 - with serious renal insufficiencies
 - in state of pregnancy

INTRODUCTION

Magneto therapy applying low frequency pulsating magnetic fields is by now recognized from numerous medical reports as a valid physical therapy, with positive biological effects, especially suitable in the treatment of bony fractures and of illness of inflammatory and chronic degenerative type.


The Extremely Low Frequency (ELF) magnetic fields seem to act on the soft structures because of micro current that are proportional to the intensity and to the variation speed of the magnetic field, being caused by the electro-motrice forces induced by the variation of the magnetic field.

To explain the interaction between the living organism and the magnetic fields one refers often to the theory of the "biological crystals" and to the crystalline-like structures of the cellular membranes, that endowed with piezo-electric propriety, would transform the applied magnetic variation in variation of their electric state and they accordingly would also determine mechanical micro modifications of the bony structures. Since the micro modifications are dependent from the direction of the induced current, the explanation of the effects reported by some patients regarding the reduction or growing acute of pains due to magneto-therapeutic treatment could be sought-after in the orientation of the magnetic field applied and consequent direction of the induced modifications.


When a therapeutic session is planned it is necessary to analyze the course of the magnetic field variation induced current: in a cylindrical solenoid the direction of the lines of the magnetic flux is parallel to the axle of the cylinder and the induced current. For instance in a femur immersed in such field, perpendicular to the above said axle and therefore perpendicular to the axle of the femur, this is true. If the same femur is immersed in a field constituted from two plain solenoids with coherent flux, the induced current would result parallel to the axle of the femur. An induced current parallel to the axle of the femur seems to give mostly positive results with reference to effect comparative experiences. That is valid for bony lesions, for instance osteoporosis or delays of consolidation, as well in presence of cutaneous substance loss, of wounds, or of soft parts.

For the same reason the stimulation of a bony lesion on the trunk or on the basin should be effected dipping the zone in a cylindrical solenoid to get induced current that are parallel to the object of treatment osseous structures.


The effect of the magnetic field mediate from the generation of electro-motrice induced forces would explain the difference of positiveness of the results reported by different researchers which used magnetic fields with different variation modality and therefore inducing different electro-motrice forces also with the same field intensity.

Absorbed femur in the cylinder. The F.E.M. is perpendicular to the axis

Femur between the plates. The F.E.M. is parallel to the axis

Tension, compression

In fact as the 50 Hz or 100 Hz sine magnetic fields varie much more slowly than the steeply variable magnetic fields produced by modern electronic methods, the current induced by the sine variation of the magnetic field is of inferior intensity and the therapeutic effect is normally blander and in any cases not remarkable.

The uncorrected use of the field generators or the intrinsic limits of the used magnetic field kind could explain the failure gotten in searches reported by some authors.

Referring to experiments effected on frogs it has been observed that the electric stimulation of stumps is determinant for their growth in innervated subjects, this allow to retain that the mechanism of new growth is managed from the nervous central system which is stimulated from the electric applied impulses directly or inductively generated from the magnetic applied field. The electric stimulation through infixed electrodes in proximity of the fractures in man has been effected with success for a long time and the obtainable results with the application of magnetic fields are similar to those gotten with the use of the aforesaid electrodes.

In many cases of pseudoarthrosis magnetotherapy seems to be the only conservative practicable treatment with probability of success and the successes obtained in the stabilization of the clinical situations or in the regression of the pathologic osteoporetic syndrome depose objectively in favour of the therapy effected through magnetic varying fields.

ELF magnetotherapy is characterized by:

- Low frequency
- Negligible electric field
- Low radiated power
- Negligible thermic effects
- High tissue penetration

For what concern the general therapeutic parameters, low frequencies and intensity relatively more elevated are normally prescribed to get a biostimulating effect, while elevated frequencies and lower intensity are preferred to achieve an analgesic and/or sedative effect. The middle frequencies and intensities are applied in the arterial circulation stimulus.

The application time vary with regard to the pathology from about ten minutes for the arthreumatic forms up to few hours for illness like pseudoarthrosis, and the number of sessions from few days for treatment of inflammatory acute processes to months for illness like pseudoarthrosis.

GENERAL INFORMATION

MT PROFESSIONAL is an ELF magnetotherapeutic unit designed for in ambulatory or at patient home professional applications, it is especially proper in the therapies that require treatments repeated for an along period and/ or long duration sessions, or therapies on traumatized or for any reason immobilized patient that therefore have difficulty to reach a physiotherapeutic ambulatory.

The unit has two square wave electric current generators, with varying frequency (2- 5- 10- 20- 30- 40- 50- 60 Hz) and varying pulse duration in order to get a predetermined and selectable percentage (20 - 40 - 60 - 80 - 99) of the intensity field related to the maximum obtainable from the connected transducer. A couple of flat transducers and a cylindrical transducer can be connected at the same time to the unit.

The couple of flat transducers and the cylinder can be activated at the same time and freely replaced with other since the unit automatically recognises them and proper own characteristics to deliver the pre-set intensity.

The magnetic field polarity is labelled on the flat transducers so to allow the use of a single transducer or add the intensities produced by two transducers with coherent fields linked.

The "north" polarity of the cylinders is situated in the side where the feeding cable is.

The intensity of the desired magnetic field is reached by modulating the actuation time of the feeding current of the solenoid which current is monitored by a warning circuit managed by an electronic microcontroller that assures its precision.

The maximum intensity of the magnetic field generated by each transducer is equal to the value of the magnetic nominal field of the same transducer; intensities lower than this value can be selected by the unit controls.

In figure 4: (1) the current and the caused field has rise time dependent from the characteristics of the solenoid.

When the value of the desired field is reached, the current and the magnetic field stops. (2) The rapid variation of the magnetic field determines the current induced in the patient.



Figure 4

The generated magnetic field is of impulsive type with pulse repetition frequency presettable by the operator, while the slope characteristics of the impulse and duration of the same is dependent both from the used transducer and the required intensity of field. When the pre-set field intensity is reached the field is quickly annulled and this rapid variation of the field induces the electro-motrice force that determines the stimulating currents in the patient.

CONTROLS

1. Mains switch
2. Section TIME:
 - a. Start Key for begin the calculation of the time
 - b. Key Up to increase the figures
 - c. Key DWN to decrease the figures
 - d. Key Stop to end the treatment
 - e. Digital Time Display
3. Section FREQUENCY:
 - a. key + to increase the working frequency
 - b. key - to decrease the working frequency
4. Output CH1 (Solenoids)
 - a. CH1 enabled light warning
5. Key MODE to enable the output channels
 - a. CH1 enabled light warning
 - b. CH2 enabled light warning
 - c. CH1 and CH2 enabled lighth warning
6. Section "INTENSITY"
 - a. Key+ to increase the field intensity
 - b. Key - decrease the field intensity
7. Output CH2 plane transducers
 - a. CH2 enabled light warning.

TECHNICAL CHARACTERISTICS

Working Frequency:	2- 5- 10- 20- 30- 40- 50- 60 Hz.
Maximum intensity:	Transducer dependent
Intensity control:	20 – 40 – 60 – 80 - 99%
Time of treatment:	1- 99 min.
Active Channels:	CH1- CH2 – CH1 + CH2
Power supply:	115 or 230 V +/- 10%- 50/ 60 Hz
Fuses:	2 x T3.15 A
Power:	250 VA
Environment Temperature:	from 10 °C up to 40 °C
Relative Moisture:	from 30% up to 75%
Atmospheric Pressure:	from 70 kPa up to 106 kPa (from 10 PSI up to 15 PSI)
Protection	IP21

SIZE

Height:	150 mm
Width:	360 mm
Depth:	250 mm
Weight:	8 Kg

STANDARD COMPOSITION

1 REF GMA50300.10	MT PROFESSIONALE electronic unit
1 REF 00100.03	Mains cord 2 m long SIEMENS-IEC

COUNTER-INDICATIONS AND COLLATERAL EFFECTS

- All the treatments have to be performed under medical control.
- The magneto-therapy use is counter-indicated for the patients:
 - pace-maker carriers or other electronic installed device carriers, and in proximity of such patient when their distance is less than one meter.
 - cardiac
 - with serious unbalances of the arterial pressure
 - with serious illness of the nervous system
 - with serious renal insufficiencies
 - in state of pregnancy
- A transitory arise of the painful symptom is possible in rheumatic patient as well the establishment of feeling of dizziness and of insomnia in applications involving the rostral segment and ureteric colic in the applications to the lumbar segment.
- When the apparatus is in operation don't approach the magnets to any type of screen (monitor, TV, etc.), this could cause the distorsion of the images and of the colours.
- The applicatory procedures are under the direct and exclusive responsibility of the doctor in charge of the case.

OPERATIONAL INSTRUCTIONS

--Inspect the package to remark possible serious damages due to the carriage. If you find damages warn the supplier and the carrier before light the unit.

--Connect the mains cable to the connector set on the back panel of the unit and, after have verified that the mains corresponds to those labelled in the unit, plug the unit to the mains net that has to be provided with earth protection or at least protected through automatic differential switch conforming the safety requirements.

--Put the magnetic field generating transducers on the treatment area and assure them by a band. One transducer alone can be used if retained sufficient.

--Connect the transducers to the appropriated connectors on frontal panel.

--Light the unit by the "POWER" mains switch.

When lighted the unit predispose the treatment time, the magnetic field repetition frequency and the field intensity percentage to the values previously pre-set before switching off the unit after the previous use. If the currently desired parameters are different from those suitable to the lighting:

-- Predetermine the desired treatment time by means of the buttons + and - in the timer section.

-- Select the desired working frequency by means of the buttons + and- in the "FREQUENCY" section.

--Select the desired magnetic field intensity percentage by means of the buttons and- in the "INTENSITY" section.

--Select the channels to activate by means of the button "MODE"

--Start the treatment by means of the "START" button in the timer section located on the left zone of the frontal panel.


--When the pre-set treatment time has elapsed, the generation of the magnetic field by the unit will automatically stop.

If during a session it wants interrupt the treatment temporarily, it needs to press the button "STOP" in the timer section; to resume the treatment press the button "START": the treatment will go on up to the pre-set time attainment noting the elapsed time before the interruption. If one desire interrupt a treatment in course definitely, press twice consecutively the button "STOP" in the timer section.

If the pre-set treatment time is 00 the treatment once on can be terminated only by hand pressing twice the button "STOP" in the timer section, this allows to effect treatments with application time exceeding 99 minutes.

TRANSDUCERS PLACEMENT

1. Linked coherent field: for in depth treatment. Polarity: North.
2. Linked not coherent field: for superficial treatments.
3. For contemporary treatments on two areas.
4. Like figure 1 with polarity: South.
5. Treatments with an only plate. Polarity: South.
6. Treatments with an only plate. Polarity: North.

MAINTENANCE

The preventive maintenance consists in the periodic cleaning that has to take place on regular base to the purpose of maintain the unit in the best working and appearance conditions.

The cleaning frequency depends on the environmental conditions in which the unit is used.

The dust could be removed with a cloth dampened with a house cleaning solution.

The malfunction problems can be of mechanical or electronic type: **MT PROFESSIONALE** incorporates some safety circuits that make practically impossible the break-up of the unit due to operational error the user.

Should a malfunction be detected it is suggested to check firstly the control setting for possible setting mistakes.

Test the mains voltage that has to power the unit, it must not exceed the label limits and check the unit to wall connection, it must be reliable.

If the unit is correctly connected to a functioning mains socket and put in the on condition does not light, check the fuses and if blow out replace with them with delayed fuses of the nominal value of 3.15A (T 3.15A). If after this substitution the unit continues to don't light or if it lights only for few minutes desist from the attempt of restoring and send the unit to a qualified center for technical assistance.

If the unit lights but the magnetic field doesn't seem to be present, when tested with a small permanent magnet, it needs to first of all check the transducers connection, if the connection is fine switch off the unit and, if resistor tester is available, disconnect the transducers from the unit and test the resistance: few ohm is the normal value: in any case less than 50 ohm. In consequence to the measurement result try the reparation or send the unit to a qualified center for technical assistance.

In any case don't open the unit cabinet since this operation if effected by personal not particularly experienced and authorized could involve from one part risks to the safety of the person that would be exposed to electric tensions which could result harmful, and from the other part the loss of the guarantee given by the manufacturing firm and/ or the supplier of the unit.

APPLICATIONS

The principal effects of the ELF magnetic therapeutic fields are:

--antiaedemic

--antalgic

--antiflogistic

--regenerative

main fields of application are:

***Traumatology :** pseudo-arthritis and delay of consolidation of the fractures, fractures with risk of pseudo-arthritis, lesions of the soft parts (bruises, distorsions, burns, etc.), district algodystrophy

***Rheumatology:** degenerative rheumatisms, inflammatory rheumatisms, fiber-mialgic rheumatisms, bony antiseptic necrosis

***Angiology:** peripheral arteriopathy, sore and ulcers

***Neurology:** muscle-tensive head-ache, post-herpetic neuralgia

***Dermatology:** psoriasis, ittirosis

***Endocrinology:** senile and post menopausal osteoporosis

Low frequency magnetic pulsating fields could be useful carriers for topic cosmetic-pharmaceutical products in order to better their absorption and therefore are suitable in aesthetical medicine as helper in the treatment of the located lipodystrophy in as they increase the vascularization and the oxygenation of the tissues.

THERAPEUTIC INDICATIVE PROTOCOLS

The applicatory procedures are under the direct and exclusive responsibility of the doctor in charge of the case.

The values of the parameters pointed out in this chapter are not critical but only indicative.

To get the actual field intensity it is necessary consider the selected intensity percentage on the unit and of the nominal intensity of the used transducer

	Frequency (Hz)	Intensity (Gauss)	Time of treatment (minutes)
Apoplexy	5	50	90
Abdominal colics	60	50	60
Acute sinusitis	50	50	30
Arthritis	30	50	30
Arthrosis	30	50	30
Asthenia	40	50	30
Bronchial Asthma	10	30	30
Bronchitis	20	50	30
Burns	30	80	60
Bursitis	20	80	90
Cellulitis	10	60	30
Cephalea	50	50	50
Chronic Sinusitis	5	40	60
Decubitus ulcers	30	80	60
Distorsions	30	80	90

	Frequency (Hz)	Intensity (Gauss)	Time of treatment (minutes)
Dyspepsias	10	80	30
Epicondylitis	40	80	60
Fractures	30	50	90
Herniated intervertebral disk	5	80	90
Infections	10	80	90
Insomnia	50	30	60
Lumbagos	50	60	90
Luxations	10	80	60
Migraine	60	20	10
Muscular atrophy	10	80	240
Muscularis pains	60	70	30
Nephritis	50	60	60
Nephrosclerosis	5	20	90
Nephrosis	10	30	90
Neuralgia of the limbs	60	60	30
Neurovegetative disorders	50	20	20
Obstructing arteriopathy	10	60	60
Osteoporosis	5	80	480
Pain ghost	5	30	90
Pseudoarthrosis	10	80	480
Rheumatic illness	60	30	60
Stomach-ache	50	50	60
Tendinitis	20	60	30
Tonsillitis	50	50	30
Traumatic lesions	10	80	240

	Frequency (Hz)	Intensity (Gauss)	Time of treatment (minutes)
Trigeminal neuralgia	50	20	20
Ulcerative colitis	30	80	60
Wounds	30	50	90

GRAPHIC SYMBOLS MEANING

Conforming to the MDD 93/42/EEC


Symbol of protection against the direct and indirect contacts: Type BF (symbol IEC 878-02-03) (CEI EN 60601-1)

Attention: consult the attached documentation (symbol IEC-ISO 348) (CEI EN 60601-1)

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