
GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403
Export: tel. +39 02 953854209/221/225 fax +39 02 95380056
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Bilancia Digitale Pegaso
Pegaso Digital Scale
Balance Numérique Pegaso
Báscula Digital Pegaso
Balança Digital Pegaso

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO

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FEATURES

The GIMA Pegaso Digital Scale is a device used to measure the height and weight of the body. It has a large digital display and an internal rechargeable battery, hence it can be used both connected or disconnected from the mains.

Maximum capacity: 200 – 250 kg

Minimum capacity: 2 kg

Precision: 100 g

Display: LCD

Dimension of weighing plate: 372 cm X 272 cm

Dimensions: 535 cm X 275 cm X 940 cm

Weight: 13 kg

Operating temperature and humidity: 5°C-40°C Max 85% RH

Power supply: 6 V Battery

Not certified product class III

PRESCRIPTIONS

Do not use the equipment in case it is damaged. Apply to your retailer.

Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

UNPACKING

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.

They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

INSTALLATION

The product is mostly factory pre-assembled. Therefore assembly operations are quite simple (see illustration):

- Fix the column to the base with the supplied screws
- Connect the connectors to the base of the column
- Insert the altimeter in the seat on the scale's column
- Fix the altimeter in the hole in the plate
- Insert the battery charger plug in the outlet on the side of the panel
- Connect the battery charger to the mains to charge the scale
- Recharge the scale every time the AC indicator lights up

OPERATION

The scale is fitted with an ON/OFF switch. Verify that the switch is set to "1" before any operation. Before the person steps onto the scale it is necessary to lift the cursor of the altimeter beyond the apparent height of the person.

After the person has stepped onto the scale, slowly lower the cursor until it touches the head.

Caution: Do not maintain the cursor in horizontal position because it may generate dangerous situations.

Function of buttons

- kg/lb** Selection of weight in kg or lb
- F** Back light function
- Z/T** Reset button to set the value to zero or set a tare value
- HOLD** by pressing this button, while weighing the person the value is blocked on the display. If the button is pressed again, the button deactivates the function. An indicator shows when the function is active
- ON/OFF** ON/OFF switch

- Switch on the scale using the ON/OFF switch, the indicator will run the initialization control and will show the figures from 000000 to 999999, at the end of the process it will stop in weighing mode. If the weight on the plate does not correspond to zero, the scale automatically activates the reset function until it reaches zero. If the scale does not run the operation automatically, it is necessary to manually adjust or recalibrate the scale.

- If the weight is lower than 2% of the total capacity the ZERO indicator will light up, press Z/T.

SETTING THE AUTOMATIC SHUT OFF FUNCTION AND ADJUSTMENT

- The scale has an automatic shut off function for energy saving purposes and switches off when it is not used for a few minutes
- Switch the scale on using the «on/off» button
- Press both «Z/T» and «F» to enter the automatic shut off function
- Use «F» to choose after how many minutes (from 1 to 9) the scale will shut off automatically

- Press «HOLD» to confirm
- The following message will appear «CAL» (scale calibration)
- Press «HOLD» to confirm or «F» to exit
- Position the maximum weight «200 kg» on the scale
- Enter 200 on the display using «F» to set the value , «kg/lb» to move to the left and «Z/T» to move to the right
- Confirm with «HOLD»
- If «ERR» appears it means that the value entered for the adjustment of the scale is too low
- If «OF» appears it means that the weight is too high and the scale is in overload

LIST OF ERROR MESSAGES

Err1	The A/D value is too low during adjustment
—OF—	Overload
O-H	The zero point is too high when the scale is switched on
O-L	The zero point is too low when the scale is switched on
P-L	Battery low

MANINTENANCE

The product does not require any special maintenance. The technical life of the scale depends on its use. Do not use water or liquid products to clean scale, juts use a soft cloth.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from.

Products sent to GIMA will be rejected.

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