GIMA

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SUPER VEGA BATTERIA is a device working with 230V ~ AC 12V DC mixed power supply.

Easily transportable from one hospital ward to another, particularly suitable for tracheotomized patients, minor surgical applications and post - operative therapy at home to be used for the nasal aspiration, oral aspiration, tracheal aspiration of the body liquids (mucus or catarrh) in the adult or in the children. Easily portable equipment designed for virtually continuous use, thanks to an electronic system for managing the electricity supply with in addition to monitoring the functioning of the internal battery, also enables the equipment to be used at a voltage of 230 V ~ even when the battery is being re-charged. A series of LED on the front control panel signals whether the machine is operating, as well as the charge and re-charge status of the batteries. Made of highly heat - resistant, electrically insulated plastic material in conformity with the latest European safety standards. Supplied with a polycarbonate autocavable jar complete with overflow valve. Equipped with aspiration regulator, located on the front panel.

GENERAL WARNING

READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

ONLY HIGHLY QUALIFIED STAFF USE RESERVED

**THE INSTRUMENT MUST NOT BE DISASSEMBLED
FOR A TECHNICAL SERVICE ALWAYS CONTACT GIMA S.p.A.**

IMPORTANT SAFETY RULES

1. On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and / or peeling of the power supply cable.
In these cases don't connect the plug to the electric socket. Carry out these controls before each use.
2. before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected.
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions.
4. Respect the safety regulations indicated for electrical appliances and particularly:
 - Only use original accessories and components;
 - The device can be used only with the bacteriological filter;
 - Never immerge the appliance into water;
 - Position the appliance on flat stable surfaces;
 - Position the device in a way that the air inlets on the back aren't obstructed;
 - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
 - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
 - Keep off the reach of children or not capable people without supervision;
 - Don't leave the appliance connected to the power supply socket when not in use;
 - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
 - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources.
5. For repairs, exclusively contact GIMA S.p.A. technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device.
6. **This medical device must be destined exclusively for the use for which it has been designed and described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force.
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents.
8. Instrument and accessories discharging must be done following current law regulations in every country of use.
9. The lead battery integrated in the device is not to be considered as an ordinary domestic waste. Such a component must be disposed of in a specific collection centre in order to be recycled.

DISPOSAL OF WASTE BATTERIES - (Directive 2006/66/EC)

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery.




The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

TECHNICAL CHARACTERISTICS

TPOLOGY (MDD 93/42/EEC)	Class IIa Medical Device
MODEL	SUPER VEGA BATTERIA
UNI EN ISO 10079-1	HIGH VACUUM / LOW FLOW
POWER FEEDING	230V~ / 50Hz or Internally powered equipment (PB Battery 12V 3Ah) or 12Vdc (1.9A) from car battery
POWER CONSUMPTION	50VA
FUSE	F 1 x 630mA 250V
MAXIMUM SUCTION PRESSURE (without jar)	-75kPa (-0.75 Bar)
MINIMUM SUCTION PRESSURE (without jar)	Less - 40kPa (-0.40 bar)
MAXIMUM SUCTION FLOW (without jar)	16 l/min
WEIGHT	4.9 Kg
SIZE	350 x 210 x 180mm
BATTERY HOLDING TIME	45 ÷ 50 minutes
BATTERY TIME CHARGE	120 ÷ 150 minutes
WORKING CONDITION	Room temperature: 5 ÷ 35°C Room humidity percentage: 30 ÷ 75% RH Altitude: 0 ÷ 2000m s.l.m.
CONSERVATION CONDITION AND TRASPORT	Room temperature: -40÷ 70°C Room humidity percentage: 10 ÷ 100% RH

SYMBOLS

	Class II isolation equipment
CE 0123	CE marking in conformity with EC directive 93/42/EEC Manufacturer: CA.MI. di Attolini Mario & C. snc Via Ugo La Malfa nr.31 – 43010 Pilastrò (PR) Italy
	Warning, consult the instruction manual
	To Preserve in place coolness and dry land
	Conservation temperature: -40 ÷ 70°C
	Type B equipment
	Fuse
	Alternate Current
	Battery
	Direct Current
Hz	Mains Frequency
I	ON
O	OFF

ACCESSORIES SUPPLIED

DESCRIPTION
COMPLETE ASPIRATION JAR 1000cc
CONICAL FITTING
TUBES SET 6 mm x 10 mm
ASPIRATION PROBE CH20
ANTIBACTERIAL FILTER

The filter is produced with (PTFE) hydrophobic material witch prevents fluids entering the pneumatic circuit. When the filter is wet, it's not possible to use the unit therefore the filter should be changed immediately. In case of possible contamination or discolouration, change the filter immediately. Don't use the suction unit without the protection filter fitted. If the suction unit is used in an emergency or in a patient where the risk of contamination is not know the filter must be changed after each use. Available under request with different versions with complete jar 2000cc.
The suction probe must be changed at the end of each application.

WARNING: Suction tubes for insertion in the human body purchased separately from the machine should comply with ISO 10993-1 standards on material biocompatibility.

Guidance and manufacturer's declaration – Electromagnetic Emissions (as request by regulation EN 60601-1-2:2001)		
The surgical aspirator SUPER VEGA BATTERIA is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator SUPER VEGA BATTERIA should assure that it's used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator SUPER VEGA BATTERIA only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The surgical aspirator SUPER VEGA BATTERIA can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes or environments feeds to you from batteries.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – Electromagnetic Emissions (as request by regulation EN 60601-1-2:2001)		
The surgical aspirator SUPER VEGA BATTERIA is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator SUPER VEGA BATTERIA should assure that it's used in such an environment.		
Immunity Test	Compliance	Electromagnetic environments - guidance
Electrostatic discharge (ESD) IEC/EN 61000-4-2	± 6kV on contact ± 8kV in air	Floors should be wood, concret or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	± 2kV power supply	Mains power quality should be that of a typical commercial environment or hospital or environments feeds to you from batteries.
Surge IEC/EN 61000-4-5	± 1kV differential mode	Mains power quality should be that of a typical commercial environment or hospital or environments feeds to you from batteries.
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	5%U _T for 0.5 cycle 40%U _T for 05 cycle 70%U _T for 25 cycle <5%U _T for 5 sec	Mains power quality should be that of a typical commercial environment or hospital If the user of the surgical aspirator SUPER VEGA BATTERIA request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Conducted Immunity IEC/EN 61000-4-6	3Vrms 150kHz to 80MHz (for appliances that aren't life - supporting)	-
Irradiated Conducted IEC/EN 61000-4-3	3V/m 80MHz to 2.5 GHz (for appliances that aren't life - equipment)	-
Note U _T is the value of the power supply voltage		

CLEANING OF ACCESSORIES

Wearing disposable latex gloves, proceed to clean the accessories as follows:

1. Grip the jar with your hand and turn the cap in a counter-clockwise direction.
2. Remove the red cage covering the float located in the jar cap (conical part with 3 openings)
3. Once the cage is removed, the float with the red seal comes out of the jar cap.
4. **To disinfect the accessories, use either denatured alcohol or hypochlorate-based solutions, easily purchasable from chemist shops.**
5. Once cleaning has been completed, insert the float in the cage, with the seal facing towards the cage opening.
6. Insert the opening of the cage into the seat in the jar cap.
7. To facilitate the cleaning of the jar, fill it 1/3 full of ordinary water. The aspirated liquid, thus diluted, will then be easier to remove.

It is possible autoclave the accessories cover and jar: insert the parts in autoclave and to effect a cycle of sterilization having care to position turned upside-down the graduated jar (with turned fund toward the tall one). After the sterilization and the cooling to the temperature the components' environment, to verify that this last don't result damaged; disassembled therefore the container for inhaled liquids following the inverse operations to the dismantlement.

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C.
The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C.
The device is ready for a new employment now.

During the sterilization operations the operator must use disposable gloves in latex, shirts and oportune mask not to enter contact with possible substances contaminating.

DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER

RULES FOR RETURNING AND REPAIRING

COMPLYING WITH THE NEW EUROPEAN RULES, GIMA S.p.A. INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE. THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.

Every returned instrument will be hygienically checked before repairing. If GIMA finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of NOT REPAIRED INSTRUMENT, accompanied by an explanation letter.

GIMA will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, GIMA will substitute the instrument, only if a SALE RECEIPT and STAMPED GUARANTEE accompany the same.

GIMA is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting.

We also request to specify the kind of fault, in order to speed up repairing procedures.

To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use.

Always specify the fault encountered so that GIMA can establish whether it falls into the category of the faults covered by the guarantee.

MAINTENANCE

The SUPER VEGA BATTERIA suction equipment does not need maintenance or lubrication.

It is necessary to check functioning and instrument before every use. Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use. Connect cable to electrical network and turn switch on.

Close the aspirator outlet with your finger and with suction regulator in maximum vacuum position check that the vacuum indicators reaches – 75 kPa (-0.75 bar) minimum (internal battery). Rotate the knob from right to left and check the aspiration regulating control.

The vacuum indicator should go down – 40 kPa (-0.40 bar). Verify that loud noises are not present, these can indicate wrong functioning.

The equipment may also be used with the power cord hooked up to the electricity mains (the yellow light switches over to the fixed green light).

Ensure that the green light comes on and that, after a few seconds it switches over to the yellow light (meaning that the charged battery is operating).

A protection fuse (**F 1 x 630mA 250V**) reachable from exterior and it situated in the plug protects the instrument.

For fuse replacing, always the type and the range.

The device is made up of a lead battery which cannot be accessed by outside. In order to replace it, consult the technical staff authorised by the manufacturer.

Fault type	Cause	Solution
1. Red light on	Battery run down	Hook up the power cord to the electricity mains, positioning the equipment power switch on 0. Charge it for 8 hours, until the yellow light switch over to the green light.
2. Red light on	Battery run down	The device may be also be supplied by an power cord hooked up to the electricity mains: the yellow light comes on and then switch over to the green light.
3. No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
4. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
5. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into it's place
6. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
7. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
8. No aspiration due to flow leakage of mucus	Filter blocked	Replace filter
Faults 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8	None of the remedies has achieved the desired results	Contact the seller or GIMA After-sales Assistance Service

If the overflow security system it's activated, don't proceede with the liquid aspiration.

If the overflow security system doesn't work there are two cases:

1° case – If the overflow security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.

2° case – If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to GIMA S.p.A. technical service.

The manufacturer will provide upon request electric diagrams, components list, descriptions, setting instructions and any other information that can help the technical assistance staff for product repair.

BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT GIMA S.p.A. TECHNICAL SERVICE.

GIMA S.p.A. DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.

INSTRUCTIONS

- Connect the short silicon tube **3**, with antibacterial filter **1**, to the suction connector **4**. The other tube, with one end connected to the antibacterial filter with the other end to jar's lid connector **2** where has been fixed the red float. When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to the avoid liquid penetration inside the device. The device must be used on a plan of horizontal operation.
- Connect the long silicon tube **5** to the other jar's lid connector **6**
- Connect the other end of the long silicon tube **5** to the probe plastic connector **8** then connect the suction probe **9** to it.
- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch **7** on position **I** to start suction
- Unscrew the jar's lid and fill the jar 1/3 full or ordinary water (this for an easy cleaning operations and an rapid reaching of the functionally vacuum) then rescrew the lid on the jar correctly.
- Once finished push switch on **0** position and unplug.
- To extract the accessories and start with cleaning.

WARNING: Ensure that the IN marker on the filter is on the side facing the collection jar lid and fitted into the "VACUUM". A wrong connection causes immediate destruction in case of contact with sucked liquids.

I prepare him it results to have the battery inside discharge, once the equipment to be used, to be able to load it is necessary to connect the cable of the device to the electric network 230V 50Hz for around 120 + 150 minutes with general interrupter of the out same device. The autonomy of the battery to full position is of around 45 + 50 minutes with operation in continuous. In the case in which instead the device is furnished with internal battery discharge, the time of recharge of the internal battery it is of around 8 hours. Then the cycle of use to 12V has to be 45 minutes ON - 120 minutes in recharge. The function of the three lights set in side position on the suction top housing is that to furnish the state of position of the battery. With battery totally position will be turned on the green light, this light it would be had to turn on as soon as we give tension to our equipment, after few seconds the yellow light mail is turned on in the middle of the three lights (intermediary functioning battery) while the red light is meaning that the battery is discharged.

The manufacturer cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.

Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC Directive and its normative.

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

At the end of its working life, the product must not be disposed of as urban waste.

It must be taken a special local authority differentiated waste collection or to a dealer providing this service.

Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant saving in energy and resources.

As a reminder of the need to dispose of separately the product is marked with a crossed-out wheeled dustbin.

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