GIMA

PULSOXIMETRO
PULSE OXIMETER
PULSOXYMÈTRE
PULSOXIMETRO

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FEATURES

Product description

COMDEK MD-600P portable pulse oximeter (cod. 34290) is an instrument with high accuracy and reliable performance.

This hand-held unit is easy to operate with the functions of adjustable alarm settings for oxygen saturation and pulse rate.

The product presents the following features:

- high accuracy
- quick assessment of Oxygen saturation and Pulse Rate
- SpO₂ and Pulse Rate display
- adjustable SpO₂ and Pulse alarm settings
- clear numerical LED display
- compact, light weight and easy to operate
- reliable quality and excellent sales service

Specifications

Oxygen Saturation Range: 0% to 100%

Pulse Rate Range: 30 bpm to 250 bpm

Display:

Oxygen saturation, Pulse Rate, Pot, Prb, Pulse search indicator & Low battery indicator

SpO₂ Accuracy: ±2% between 100% to 80%
±3% between 70% to 65%
below 65%, not specified

Pulse Accuracy: ±1% of full scale

Default Alarm Settings:

- SpO₂ Upper Limit:: 100%
- SpO₂ Lower Limit : 85%
- Pulse Rate Upper Limit: 150 bpm
- Pulse Rate Lower Limit: 30 bpm

Power Requirements:

- DC 6 Volts by 4 pieces AA size batteries
- AC AC adapter output DC 6 V, 600 mA
- rechargeables AA batteries (optional)

Operating Temperature:

-20°C to 55°C

SpO₂ Sensors:

Finger clip sensor & wrap probe (optional)

Dimensions:

172 mm x 90 mm x 36 mm

Weight:

240g

Field of application

This is a pulse oximeter for the 21st century: reliable, portable and affordable. It is designed to the fields of hospitals, clinical surgery and home care applications.

PRESCRIPTIONS

Accident prevention

The product shall be used by qualified personnel only.

Do not use the equipment in case it is damaged. Apply to your retailer.

Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

TRANSPORT

Transport procedures

It is not necessary to follow particular transport procedures. It is however recommended to carefully handle it and avoid shocks. The equipment is delivered in an imitation leather case allowing a comfortable shoulder transport.

Unpacking

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.

They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

INSTALLATION

Installation and connections

The pulse oximeter is delivered with the following components: a practical imitation leather case, 4 “AA” alkaline batteries, AC adapter with 6 V output, 600 mA, a finger sensor for heart beat survey. To feed the equipment with 4 “AA” batteries, open the apposite compartment on the back and insert the batteries. Pay attention not to invert the direction of poles (see the picture inside the battery compartment) and then replace the cover.

To feed the equipment with an AC adapter, insert the adapter jack into the corresponding socket on the equipment base and plug in the adapter.

With the exception of rechargeable batteries, remove common or alkaline batteries when using the AC adapter for charging.

Insert the finger sensor connector into the pulse oximeter (when the sensor is not inserted, the display views the following writing: **Prb**). It is sufficient to make the grooves present on the connector and on the socket match.

OPERATION

Instrument description and symbols

Select Pushbutton for the selection of SpO 2 Max and Min values and PULSE/min.

Pushbutton to increase SpO 2 Max and Min values and PULSE/min.

Pushbutton to decrease SpO 2 Max and Min values and PULSE/min.

Turn on/off audio signals

Indicator for disabled audio signals.

Run-down battery indicator. When the battery charge is low, the battery indicator lights up and gives out two beeps. In this case it is necessary to replace the batteries. At this point the run-down battery indicator turns off.

Prb Sensor not connected with the equipment

Pot No patient connected

Lo Low signal

Operating instructions

To turn the equipment on, press the green ON/OFF pushbutton. To confirm power-up, an acoustic signal is released and **PWR** pilot light turned on. The right display views the following writing: **Pot** (no patient connected).

Insert the finger of the patient into the examination sensor. It is possible to use any finger, however it is advisable to select the finger suitable to the sensor dimensions. (The PULSE OXIMETER can be fitted with different sensors, follow the instructions).

In order to perform an accurate detection of the value, the patient shall keep still.

Wait 6 seconds and the sensor starts up measuring. The left display shows the value of blood oxygen saturation rate (SpO 2 %) and the right display shows the heart rate (Pulse/min).

COMDEK PULSE OXIMETER is provided with a sound signal warning that the values detected by the sensor are higher or lower than the ones stored in the equipment.

Default alarm settings are:

SpO 2	Max	100%
SpO 2	Min	50%
Heart rate	Max	250 bpm
Heart rate	Min	30 bpm

To change these values, proceed as follow:

- **SpO 2 alarm max. value**

Press **Select** once, the red led (SpO 2 Min) lights up and a value blinks on the left display. This value represents the **maximum** value for oxygen saturation (SpO 2).

To change this value, press key ▲ to increase it and key ▼ to reduce it.

- **SpO 2 alarm min. value**

Press **Select** two times, the red led (SpO 2 Min) lights up and a value blinks on the right display. This value represents the **minimum** value for oxygen saturation (SpO 2).

To change this value, press key ▲ to increase it and key ▼ to reduce it.

- **Heart rate alarm max. value**

Press **Select** three times, the red led (Pulse Max) lights up and a value blinks on the left display. This value represents the **maximum** value for heart rate (Pulse/min). To change this value, press key ▲ to increase it and key ▼ to reduce it.

- **Heart rate alarm min. value**

Press **Select** four times, the red led (Pulse Max) lights up and a value blinks on the right display. This value represents the **minimum** value for heart rate (Pulse/min). To change this value, press key ▲ to increase it and key ▼ to reduce it.

When it stops blinking, the equipment returns to the original conditions.

The new alarm settings will be automatically saved.

When the batteries run down or when the equipment is disconnected from power supply, the original alarm settings are restored.

During measurement, the pulse oximeter display shows an intermittent scale of 10 leds indicating the search for beating. As seconds go by, the scale grows longer and becomes stable. When the signal is good, the search bar should show at least 4 leds.

If the loud speaker is on, the equipment gives out intermittent audio signals. Press  pushbutton to stop the sound and the red led  on the display lights up. Press the pushbutton again and the sound starts up again.

The alarm keeps on sounding for a minute and then it stops.

Switching off

After the measurement has been taken, turn the equipment off by means of ON/OFF pushbutton. A beep is given out and then the equipment stops.

MAINTENANCE

General maintenance instructions

The product has been conceived to last forever; it is resistant to most part of chemicals except for methyl ethyl ketones (M.E.K.), diluents and nail polish removers.

To clean the instrument, make use of a cloth dampened with disinfectant.

GUARANTEE

General guarantee conditions

Thank you for having acquired our product.

This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, doppler probes, SpO₂ sensors, lamps, batteries, electrodes, handpieces, washers and packings, resistances, etc.

No compensation will be given for loss of use of the product.

Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.

The guarantee will be invalidated if the registration number has been removed, cancelled or altered. Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

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DET NORSKE VERITAS

FULL PRODUCT QUALITY MANAGEMENT CERTIFICATE - EC

Certificate No. 52177-2009-CE-RGC-NA Rev. 1.0

This Certificate consists of 3 pages

This is to certify that the Quality Management System of

Comdek Industrial Corp.

Taipei, Taiwan

for production and final product inspection/testing of

Pulse Oximeter / Patient Monitor

has been assessed with respect to

the conformity assessment procedure described in Article 11.3.a and Annex II excluding section 4 (Module H) of Council Directive 93/42/EEC on Medical Devices, as amended, and found to comply

Further details are given overleaf

Place and date:

Høvik, 11 May 2009

This Certificate is valid until:

16 February 2014

For DET NORSKE VERITAS CERTIFICATION AS
Norway

Marianne Spøren
Certification Manager

Notified Body No.:
0434

Aud Løken Eiklid
Technical Reviewer

Notice: The certificate is subject to terms and conditions overleaf. Any significant changes in design or construction may render this certificate invalid.

If any person suffers loss or damage which is proved to have been caused by any negligent act or omission of Det Norske Veritas, then Det Norske Veritas shall pay compensation to such person for his proved direct loss or damage. However, the compensation shall not exceed an amount equal to the amount that the client paid for the service in question, provided that the maximum compensation shall never exceed USD 100,000. In no event shall Det Norske Veritas be liable for any consequential, indirect, or other damages, including lost profits, arising out of or from the use of the services of Det Norske Veritas.

Cert. No.: 52177-2009-CE-RGC-NA
Rev. No.: 1.0
Project No.: PIUC-13290-2007-PRC-RGC

Jurisdiction

Application of Council Directive 93/42/EEC of 14 June 1993, adopted as 'Forskrift for Medisinsk Utstyr' by the Norwegian Ministry of Health and Care Services.

Certificate history

Revision	Description	Issue Date
	Original certificate	2004-02-16
1.0	Recertification	2009-02-16

Products covered by this Certificate

Product Description	Product Name	Class
Pulse Oximeter	- Portable Pulse Oximeter MD-600P/610P/MD-620P/MD-630P/MD-650P/ MD-660P/MD-670P/MD-680P - Oxi-Capnography MD-660P - Desktop Pulse Oximeter MD-700/MD-710/MD-720/MD-730/MD-750/MD-760/ MD-770/MD-780	IIB
Patient Monitor	MD-800/MD-810/MD-820/MD-830/MD-850/MD-860/ MD-870/MD-880	IIB
Accessory	- Sensor: D001/D001-A/D002/D003/R002/R003/R004/R004-i/ R005/RA001/R008/E009 - Probe: HD-T2252-AG/HD-T2252-AS	IIB

The complete list of devices is filed with the Notified Body.

Sites covered by this certificate

Head office: 9F-1, No. 3, Yuan Qu Street, Nan-Kang, Taipei 11503, Taiwan, R.O.C.

Manufacturer: 7F, No. 1, Alley 6, Lane 235, Pao Chiao Rd., Hsin Tien, Taipei, Taiwan R.O.C.

EU Representative

Comdek Computer System, Schwendstrasse 32, CH-8486 Rikon, Switzerland