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Instruction Manual for Pulse Oximeter (REF LTD805 - LTD807 pediatric)


1. General Description

Haemoglobin Saturation is percentage of HbO₂ (Oxyhemoglobin), compounded with oxygen, in haemoglobin (Hb). In other words, it indicates the saturation of HbO₂ in blood. It is a crucial physiological parameter for Respiration System. Many respiratory diseases can lead to decrease of SPO₂ saturation. Moreover, such factors as organic automatic adjustment malfunction derived from anesthesia, side effect of major operation and injures aroused from physical check can also cause problems about SPO₂ saturation variation, which contributed to such adverse symptoms such as faint, vomit and feel weakness and so on. What's more, if can not get the right remedy, it can lead to dangers to human's life. The Fingertip Pulse Oximeter with features of small dimension, low power consumption, easy-to-use and convenience for portable. As long as you put your finger in the device which read data through a light-sensitive probe, then the associated reading will be displayed on the indicator. It is proved by clinical trial Pulse Oximeter have outstanding precise and could repeat to measure the SPO₂ more accurately and stably.

2. Measurement principle

The principle of the pulse oximeter is based on the fact that the colour of blood depends on the quantity of the Oxyhemoglobin. While the Oxyhemoglobin is red, the colour of the non-saturated Hemoglobin tends to blue. The hemoglobin changes its structural configuration when it gets involved in a chemical reaction; each of these configurations presents a different type of light reflection. The wave length of 660nm corresponds to red colour of the Oxyhemoglobin (HbO₂), the non-saturated Hemoglobin corresponds to the wave length of 940nm. Both of the wave lengths pass through a blood sample and show the percentage of HbO₂. This percentage corresponds to the value of saturation (SpO₂).

Diagram of Operation Principle

1. Glow and Infrared-ray Emission Tube
2. Glow and Infrared-ray Receipt Tube

3. Technical Parameters

- 3.1 Display Format: LED display
Hemoglobin saturation Display SPO₂: 0-99%
Pulse Ratio Display: 30-254 BPM
Pulse rate display: barograph
Battery Voltage Low Indication:
- 3.2 Power Requirements: 2 AAA size 1.5V 3.3
- 3.3 Power Consumption: less than 40mA
- 3.4 Measurement Accuracy:
SPO₂: +2% between 80%-99%,
+3% between 70% - 80%, <70% -no definition
Pulse rate: +2 BPM between 30-235 BPM
- 3.5 Dimension:
LTD805 length: 58mm LTD807 length: 50mm
width: 30mm width: 29/26mm
height: 34mm height: 28mm
weight: 50g weight: 50g
(including two AAA batteries) (including two AAA batteries)
4. Product Properties
- 4.1 Easy to use.
- 4.2 Small dimension, light (total weight is about 50g including batteries) and portable.
- 4.3 Low power consumption; included two AAA batteries can last up to 30 hours.
- 4.4 Low voltage warning will be indicated on display when battery voltage is too low and normal application of the Oximeter might be influenced.
- 4.5 No signal can be tested, device will power off automatically in 8 seconds.

5. Application field


You can use it in the field such as family, hospital, and sports (do not use the device during exercises, but can use it before or after exercises).

Attention: The product is not suitable to monitor patient continuously.

6. Instructions

- 6.1 Installing two AAA batteries into battery cassette before closing its cover.
- 6.2 Open the clamp shown as in the picture below:
- 6.3 Put one of your fingers into rubber hole of the Oximeter (it is better to let your finger touch the bottom) before releasing the clamp.
- 6.4 Press the switch button for one time on front panel.
- 6.5 Do not move your body at the same time.
- 6.6 Read correspondent data from display screen.


6.7 Six display modes. After turn on the oximeter, each time you press the power switch, the oximeter will change display mode, there are 6 display modes shown as follows:


When you press the power switch for a long time (more than one second), the brightness of the oximeter will be changed by degrees, there are 10 levels on brightness; the default level is level four.

Attention: When you put finger into the Oximeter, your nail surface must be upward.


Attention: Use the medical alcohol to clean the rubber touching the finger inside of Oximeter, before and after operation. (The rubber inside of the Oximeter is medical rubber, which has no toxin, and no harmful to the skin of human being).

7. Brief Description of Front Panel

- 7.1:
 - a) Power switch
 - b) PR Bargraph
 - c) Hemoglobin Saturation (SpO₂)
 - d) Heart rate (PR)

- 7.2:
 - a) SpO₂ waveform
 - b) Low power indicator


The heart rate bargraph display corresponds with pulse rate.

8. Pieces of Accessories

- 8.1 One hanging string
- 8.2 Two batteries
- 8.3 One User manual

9. Battery Installations

- 9.1 Push the two AAA batteries into battery cassette in correct polarities.
- 9.2 Push the battery cover horizontally along the arrow shown as below.

Notes: Battery polarities must be correctly located. Otherwise, might cause serious damage to the device. Please put or remove batteries in right order, or is likely to damage the device bracket. Please remove the battery if the Oximeter will not be used for a long time.

10. Hang string Installations

- 10.1 Thread thinner end of the hang lace through the hanging hole. Thread thicker end of the lace through the threaded end before pulling it tightly

11. Maintenance and Storage

- 11.1 Replace the batteries timely when low voltage indicator is on.
- 11.2 Clean surface of the fingertip Oximeter before it is used in diagnosis for patients

11.3 Remove the batteries inside the battery cassette if the Oximeter will not be operated for a long time.

11.4 It is better to preserve the product in a place where ambient temperatures range from -10° e i 40°C and humidity range from 10 e 80%.

It is recommended that the product be kept in a dry place. A damp ambient might affect its lifetime and even might damage the product.

Notes: Please follow the articles of the local government to waste with used batteries

- 11.5 Environment requirements:
Operation Temperature: 5 - 40°C
Storage Temperature: -10 - 40°C
Humidity: 10 - 80%

12. Possible Problems and resolutions

Problems	Possible reason	Solution
Spo ₂ or PR can not be shown normally	1.Finger is not plugged correctly 2.Patient's Oxyhemoglobin value is too low to be measured.	1.Retry putting in the finger 2. Try some more times, if you can make sure about no problem existing in the product.
Spo ₂ or PR is shown unstably	1.Finger might not be plugged deep enough 2.Finger is trembling or patient's body is in movement status	1.Retry putting in the finger 2.Try not to move
The Oximeter can not be powered on	1.Powerof batteries might be inadequate or not be there at all 2.Batteries might be installed incorrectly 3.The Oximeter might be damaged	1.Please replace batteries 2.Please reinstall the batteries 3.Please contact with local customer service centre
Indication lamps are suddenly off	1.The product is automatically powered off when no signal is detected longer than 8 seconds 2. Low battery	1.Normal 2.Replace the batteries
"Error 3", "Error 4" or "Error 7"	1 Low power 2 Receiving or emission tubes shielded or damaged. 3 Current control circuit malfunctions.	1 Please change battery 2 Please contact with local customer service center

Declaration: EMC of this product comply with IEC60601-1-2 standard.

13. Signal Definitions

Signals	Definitions
	The equipment is BF type
	Attention: Please reference instruction manual before operation
%SpO ₂	Hemoglobin saturation
BPM	Pulse rate or Heart rate (Beat per Minute)
	Low Voltage Indication
SN	Serial No

14. GUARANTEE

This device has 2-year warranty from the purchase date which is certified by the receipt or invoice. The distributor reserves the right either to repair or replace necessary parts to correct defects in the materials or workmanship without any charge. The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alteration or disassembly by unauthorized individuals. In the case you need to get your device repaired and you fulfill the conditions of warranty, send the device to the nearest retailer together with the purchase receipt.

Distributed by:

MORETTI <http://www.morettispa.com>

Beijing Choice Electronic Technology Co., Ltd
Shanghai International Holding Corp.
GmbH (Europe) Eiffestrasse 80 20537 Hamburg Germany