ECG PALMARE PM10

PM10 PALM ECG MONITOR

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MONITOR ECG PORTÁTIL PM10

ECG PALMAR PM10

PM10 PALM ECG MONITOR

Manuale d'uso - User manual - Manuel de l'utilisateur
 Guia de Uso - Guia para utilização - Gebrauchsanweisung

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.
ATTENTION: The operators must carefully read and completely understand the present manual before using the product.
AVIS: Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.
ATENCIÓN: Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.
ATENÇÃO: Os operadores devem ler e entender completamente este manual antes de usar o produto
ACHTUNG: Die Bediener müssen vorher dieses Handbuch gelesen und verstanden haben, bevor sie das Produkt benutzen.

REF 33246 / PM10

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ENGLISH

Foreword

Thankyou very much for purchasing the Portable ECG Monitor.

This user manual introduces detail product information about its character, requirement, structure, performance, specification, appropriate methods of transportation, installation, usage, operation, repair, maintenance and storage, and safety measures of how to protect the operator and product. Please read details in the following chapters.

Please read the user manual carefully before using the product and strictly follow its regulations to operate. The user manual indicates the operations that users need to pay much attention to, that may lead to abnormality, or may danger to the device or human body during using. Our company will not response the security, reliability and performance for any abnormality or device and human body damage caused by not following this user manual to use, maintain and store, not provide free service for any situations above. We apologize for the ccontent in the manual is subject to change according to product upgrades without notice.

The product is reusable as a medical instrument.

Warning:

- The product is not a examination device apllied in clinical medicine, and its results can not serve as the basis for diagnosis, but cn be used as a reference for patient to take further medical treatment and reference for doctor to diagnose.
- The reliability depends on whether users are following the operation and maintenance in the user manual or not.
- All servicing and future upgrade to the device must be carried out by personnel trained and authorized by our company, and using the original fitting for maintenance.

This iser manual contains proprietary information, wich is protected by copyright. All rights reserved.

Reproduction, adaption or translation, for any part of the manual without prior written permission, is prohibited.

Our company takes the responsibilities as follows:

1. To provide qualifield products according to enterprise standard for users;
2. To provide services of installation, debugging and training according to the contract;
3. To provide one year warranty and product maintenance after warranty period according to the contract;
4. To respond user's requests in time.

Chapter 1 Notice

1.1 Generic notice

- 1) Do not use the device in locations subject to high temperatures or humidity. Use in the temperature within 5 to 40°C and humidity within 25% to 80% RH.
- 2) Do not wash the device with water.
- 3) Do not use or store the device in the following ambient conditions:
 - Near fires or open flames
 - Locations exposed to strong vibration
 - Locations exposed to strong electromagnetic fields
- 4) Do not sterilize the device in autoclave or gas sterilizer.
- 5) The device service life is 3 years. Do not throw away the device and accessories when they can't work. If the device needs to dispose, it should meet the local laws and regulations requirement.

1.2 Measurement notice

- If your skin is dry, wipe them with disinfectant alcohol or electric salve to strengthen the electric capability
- You are better to confortably sit, draw yourself up, begin to measure when the waveform level off.
- When measuring, the finger and chest electrodes should touch your skin exactly, roundly and well.

1.3 Safety notice

- 1) No sampling in the battery-charging.
- 2) Lay the device in shady and cool environment when you are not going to use it for a long period of time, and electrify per three months.
- 3) Do not use the device in the environment placed inflammable objects, such as anesthetic.

1.4 EMC notice

Please note the effect from EMC when using the device, because it can ben influenced by portable or movable high electromagnetic compatibility RF devices.

Chapter 2 Introduction

The portable ECG monitor is designed for family and individual users. It is a good helper for family members to prevent from cardiovascular disease, as it can monitor patients ECG anytime at anyplace with easy operation. The device can record, analize and display user's ECG waveform, capture the pathological ECG waveform when user happen to heart attack or other unpleasant symptoms. The ECG monitor can be conducted not limited in the hospital, which saves money from the physical check-up for users. After connected with a computer, users can print their ECG waveform, which provides data reference for doctors.

2.1 Characteristics

- 1) Handsome shape, handy operation, convenient tote.
- 2) Monitor and record real-time ECG waveform and HR anytime and anywhere.
- 3) Built-in large capability rechargeable lithium battery, continuously sample 200 ECG waveform after charged once.

2.2 Application

- 1) Occaasion: family, medical clinic and hospital. The device can't be used as a general electrocardiogram for clinical examination.
- 2) Object: people under high pressure and workload for long time, heart disease patients, middle aged and aged poeple, sub-health people.
- 3) Purpose: the device is only used for ECG monitoring and data storage. It is not a therapy equipment. Operation method is simple and less requirement for the operating personnel.

Chapter 3 Primary Technical Orders

3.1 normal work environment

- 1) Operation environment
 - Temperature: +5°C~+40°C
 - Relative humidity: 25%~80%
 - Power supply: built-in rechargeable lithium battery, voltage: 3.7V
- 2) Transportation and storage environment
 - Temperature: -40°C~+55°C
 - Relative humidity: ≤95%

3.2 Basic parameters

- 1) Calibration voltage: 1mV±5%.
- 2) Standard sensitivity: 10mm/mV±5%.
- 3) Amplitude frequency characteristic: standard: 10Hz; 1Hz-20Hz; (+0.4dB, -3dB).
- 4) Noise level: <30μV.
- 5) CMRR: >60dB.
- 6) Scanning speed: 25mm/s+5%.
- 7) Sampling rate: 250 dots/s.
- 8) HR measurement range: 30bpm-300bpm, error: ±1bpm or 1%.
- 9) Type of protection against electric shock: Internal power device.