PIEGHEVOLE

RIGIDO

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USER MANUAL
SCANGLISS 2 PLUS

INTRODUCTION

Thank you for choosing Demarta Virginio. The products have been designed and built to meet all your needs for a practical, correct and safe use. This manual contains some small suggestions for a correct use of the device of your choice and some precious advice for your safety. We recommend that you carefully read the entirety of this manual before using the item you purchased. If in doubt, please contact your dealer, who will be able to help and advise you correctly.

The manufacturer reserves the right to make changes to the device and to this manual without prior notice in order to improve its characteristics.

GENERAL WARNINGS

- Use the product in the manner indicated in this manual. If doubts persist, contact DEMARTA VIRGINIO SAS.
- Always evaluate the compatibility of the product with the characteristics of the subject together with the doctor or therapist.
- Always check the integrity of the product before each use, even when new. Do not use if damaged or worn.
- DEMARTA VIRGINIO SAS declines any responsibility for damage to property or persons resulting from improper use, not agreed with a professional figure, not corresponding to what is indicated in the user manual.
- Any accidents, anomalies, breakages must be reported to DEMARTA VIRGINIO SAS in order to monitor the quality and criticalities of the product.
- Patient handling operations must only be carried out by qualified staff.

DECLARATION OF CONFORMITY

DEMARTA VIRGINIO SAS declares under its sole responsibility that the product made and traded by DEMARTA VIRGINIO SAS comply with the applicable provisions of the 2017/745 Regulation on Medical Devices of 5 April 2017.

For this purpose, DEMARTA VIRGINIO SAS guarantees and declares under its sole responsibility what follows:

- The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
- The devices ARE NOT MEASURING INSTRUMENTS.
- The devices ARE NOT MADE FOR CLINICAL TESTS.
- The devices are packed in NON-STERILE BOX.
- The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation.

INTENDED USE

The Transfer Systems are intended for facilitating the manual handling operations of patients through thrust, traction, sliding maneuvers. They have been designed for use in hospitals, nursing homes and healthcare facilities and for home care.

DESCRIPTION

The transfer board SCANGLISS 2 PLUS allows one or more operators to transfer a patient from one surface to another with a difference in height up to 15 cm and in distance up to 20 cm. The transfer board is equipped with two fireproof ABS handles for transport. The transport handles can be removed unscrewing the fixing screws. It is radiolucent.

Maximum Capacity: 300 Kg.

There are two types available:

- Foldable.
- Rigid.

Product	Code	Dimensions (cm)	Weight (Kg)
Foldable Maxi	03.04.9000	180 x 50	3,90
Foldable Midi 150	03.04.9004	150 x 50	3,00
Foldable Midi 125	03.04.9002	125 x 50	2,70
Product	Code	Dimensions (cm)	Weight (Kg)
Rigid Maxi	03.04.9001	180 x 50	3,90
Rigid Midi 150	03.04.9005	150 x 50	3,00
Rigid Midi 125	03.04.9003	125 x 50	2,70
Mini	03.04.9010	90 x 50	1,50



Materials	
Handles	ABS Plastic
Base	Etaflon
Base Coating	Polyester
Ultra Low Friction Sheet	Nylon RIP-STOP
Sheet Maintenance and Cleaning	
	Water Temperature not Above 60°C or 140°F
	Do not Iron
	Do not Tumble Dry
	Do not Dryclean
Frame Maintenance and Cleaning	
Disinfection is possible with 5% diluted disinfectants even with alcohol base	
Storage and Disposal	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	

Removable handles in fireproof ABS for a more comfortable transport.

The transport handles can be removed unscrewing the fixing screws.

USE

- * Slightly lift the patient sideways.
- * Partially insert the article under the patient.
- * Push, with due precautions, the patient's body: the thrust provided causes the sliding of the RIP-STOP Nylon sheet and, consequently, also the lateral movement of the patient himself.

IMPORTANT! Maximum height difference up to 15 cm and a maximum distance up to 20 cm.

GENERAL DISPOSAL CONDITIONS

When disposing of the device, never use housing solid waste collection systems.

It is necessary to dispose of the device through the common ecological islands common for the recycling of used materials.

WARRANTY

DEMARTA products are guaranteed from material or manufacturing faults for 2 years from the purchasing date.

DEMARTA does not guarantee for damages caused by the following condition: natural disaster, not authorized maintenance and failure to adhere to the guidelines and instructions for use.