# MIO-IONOTENS

---

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**I.A.C.E.R. Srl**  
**Via S. Pertini 24/A – 30030**  
**Martellago (VE) ITALY**  
**Tel. 041.5401356 • Fax 041.5402684**  
**e-mail: [iacer@iacer.it](mailto:iacer@iacer.it) <http://www.iacer.ve.it>**

MNPG19 Rev. 2 del 21/03/10

Elettroterapia modello

# MIO-IONOTENS

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Electrotherapy model

# MIO-IONOTENS

<b>Summary</b>	<b>21</b>
Manufacturer	22
Declaration of Conformity	22
Specifications	22
Purpose	22
Technical specifications	23
Labelling	23
Label details	24
Symbols	24
Kit contents	24
<b>How to use</b>	<b>25</b>
Warnings	25
Electromagnetic interference	25
Side effects	25
Instructions	26
TENS and ionophoresis	27
Stored programs	29
<b>Maintenance</b>	<b>32</b>
Replacing the batteries	32
Cleaning	33
Carriage and storage	33
Disposal	33
Troubleshooting	33
Repairing	33
Spare parts	33
Warranty	34

## Manufacturer

**I.A.C.E.R. S.r.l.**

Via S. Pertini, 24/a • 30030 Martellago (VE) ITALY  
Tel. 041.5401356 • Fax 041.5402684

IACER S.r.l. is an Italian medical devices manufacturer (CE medical certificate n° MED 24021).

## Declaration of Conformity

IACER S.r.l., headquartered in Italy, via S. Pertini 24/A 30030 Martellago (VE), declares on its own responsibility that MIO-IONOTENS is manufactured in conformity with Council Directive 93/42/EEC (MDD) dated 14 June 1993 (D. Lgs. 46/97 dated 24 February 1997 “Attuazione della Direttiva 93/42/CEE concernente i dispositivi medici), Annex II as modified by Directive 2007/47/CE dated 5 September 2007 (D. Lgs 37/2010 dated 25 January 2010).

Notified Body: Cermet, Via di Cadriano 23 – 40057 Cadriano di Granarolo (BO) ITALY

MIO-IONOTENS is a Class IIa equipment, with reference to Directive 93/42/EEC (MDD), annexed IX rule 9 (and following modifications).

Martellago, 21/03/2010

Legal representative  
Mario Caprara

## Specifications

MIO-IONOTENS has the following specifications:

- Class IIa equipment (Directive 93/42/EEC, annexed IX rule 9 and following modifications);
- Class II, applied part type BF (Classif. EN 60601-1);
- Equipment not protected against liquids penetration;
- Equipment and accessories not subjected to sterilization;
- Use of the equipment is prohibited close to flammable substances or in environments with high concentrations of oxygen;
- Continuous operating mode equipment;
- Equipment not suited to be used in external.

## Purpose

Clinical purpose: Therapeutic

MIO-IONOTENS is recommended both for clinical and for domestic use.

MIO-IONOTENS pathologies applications:

- wrist articulation
- hand articulation
- shoulder articulation
- foot articulation
- ankle articulation
- knee articulation
- skeletal motor apparatus
- atrophies and muscular dystrophy
- bruises
- neuralgias
- benign lesions and muscular tears

- tendonitis

Thanks to its TENS protocols, MIO-IONOTENS is particularly suitable for pain therapy. TENS impulses can reduce and, in many patients, eliminate the pain generated from pathologies up above described.

MIO-IONOTENS has also specific ionophoresis protocols. Ionophoresis is an electrotherapeutic technique that uses continuous current to introduce drugs on pain or contracture area. The current promotes the migration of the drug ions: the drug passes through the pain area releasing its specific action. Ionophoresis has two great advantages: it avoids the administration of drugs by mouth and it treats directly the pain areas.


Ionophoresis is also used for the treatment of diseases affecting urogenital male apparatus, like IPP (Induratio Penis Plastic) or La Peyronie disease. Consult a specialist before start the therapy. Contact the manufacturer for other informations.

## Technical specifications

Power supply	4 x 1,5V AAA Size Alkaline batteries
Max consumption.	1,08W during treatment
Insulation class (CEI EN 60601-1)	II
Applied part (CEI EN 60601-1)	BF
Dimensions (mm)	118x62x25
Max output current	40mA peak-peak, 1 K $\Omega$ load for P01, P02, P03, P15, P16 100mA peak-peak, 1 K $\Omega$ load for other programs
Waveform	Biphasic compensated square wave (Monophasic square wave for ionophoresis programs P01-P02-P03)
Waveform frequency (Hz)	0,5÷1200 (adjustable 0,5÷150 Hz for user programs)
Impulse width ( $\mu$ s)	30÷300
Timer	20 or 30 min. for stored programs, 1÷60 min. or continuous for user programs


WARNING. The equipment delivers current in excess of 10mA.

## Labelling

## Label details

### Label 1

### Label 2

### Label 3

### Label 4

### Label 5

## ON OFF

## Symbols

	Warning. Read annex documents
	Product subject to WEEE regulations concerning separate waste collection of electronic equipment
	Class II equipment
	Applied part type BF
	Compliance with Directive 93/42/EEC (MDD) (and following modifications)
	Manufacturing date

## Kit contents

- n° 1 MIO-IONOTENS device;
- n° 4 1,5V AAA Size Alkaline batteries;
- n° 2 electrostimulation cables;
- n° 1 set of 4 48x48mm pre-gelled electrodes;
- n° 1 set of 4 50x90mm pre-gelled electrodes;
- n° 1 ionophoresis kit (1 elastic bend, 2 rubber electrodes, 2 sponges)
- n° 1 belt clip;
- n° 1 transport bag;
- n° 1 user manual.

## Warnings

- Take care of position and meaning of the labels on MIO-IONOTENS;
- Do not damage connection cables and avoid to roll up the cable around the device;
- Use of the device is prohibited in the presence of signs of deterioration of the device itself: read “Assistance” paragraph for details;
- Avoid the use of MIO-IONOTENS to people not educated through the reading of the manual;
- Do not use MIO-IONOTENS in damp environments;
- Do not wear metallic objects during therapy;
- It is forbidden to position the electrodes in such a way that the current crosses the heart area (e.g. a black electrode on the chest and a red electrode on the shoulder blade);
- Use of the device is prohibited with electrodes positioned on or close to injuries or cuts;
- The electrodes must not be positioned on the carotid sinuses (carotid) or genitals;
- The electrodes must not be positioned close to the eyes; make sure that the current delivered does not cross the eyeball (one electrode diametrically opposite to the other in relation to the eye); keep a distance of at least 3 cm. from the eyeball;
- Insufficiently sized electrode sections can cause skin reactions or burns;
- Do not use electrodes when damaged, even if they stick to the skin well;
- Electrodes must not be used when they no longer stick to the skin. Repeated use of the same electrodes can compromise the safety of the stimulation, in fact it can cause skin redness that can last for many hours after stimulation.

WARNING. The manufacturer is considered responsible of the performances, reliability and safety of the instrument only if:

- possible additions, changes and/or reparations are effected by authorized personnel;
- the electric plant of the environment in which MIO-IONOTENS is inserted is conforming to the national laws;
- the instrument is employed in conformity to the instructions contained in this manual.

## Electromagnetic interference

MIO-IONOTENS doesn't produce and doesn't receive interferences from other equipments. However, it's opportune to use the instrument keeping it at least at 3 meters from televisions, monitors, cellular telephones or any other electronic equipment.

## Side effects

Patient in pregnancy, tuberculosis, juvenile diabetes, viral (in acute phase) illnesses, mycosis, dermatitis, cardiopathic subjects, serious arrhythmias or pacemaker carriers, children, metallic prosthesis carriers, acute infections, epileptics (different medical prescriptions excepted).

No significant side effects are known of. Some particularly sensitive people have reported skin redness in the area where the electrodes were positioned: the redness usually disappears a few minutes after the end of the treatment. Should the redness persist please consult a doctor.

In rare cases, stimulation carried out in the evening can cause some people to experience difficulty in falling asleep. If this occurs, suspend the treatment and consult a doctor.


## Instructions

MIO-IONOTENS is a portable device for TENS and ionophoresis treatment designed for everyday use for the most common forms of muscular pain. MIO-IONOTENS has two independent channel adjustable by the user.


MIO-IONOTENS has 16 stored programs (5 ionophoresis and 11 TENS programs) and 4 adjustable programs: the user can adjust his own program, manually setting the treatment parameters.

MIO-IONOTENS has a special backlight function:

- In the normal state has the green color backlight.
- If the electrodes is in the wrong connection way it has red color backlight;
- If the output intensity is reaching 2/3 max value it has the blue color backlight.

1. MIO-IONOTENS
2. ON/OFF switch
3. Channel 1 and 2 output
4. Channel 1 and 2 intensity control
5. LCD Display
6. Programs selector
7. Set parameters
8. Increase value
9. Decrease value
10. Battery case
11. Belt clip

1. Intensity values
2. User programs
3. Stored programs
4. Output waveform
5. Low battery icon
6. Setting value
7. Key lock

## SETTINGS

### 1. ELECTRODES AND CABLES CONNECTION

Place the electrodes on the skin near the painful are (see next chapter for details), connect the electrodes to the cables and then connect the cables to CH1 and CH2 outputs on the upper side of MIO-IONOTENS;

### 2. SWITCHING ON

Switch on MIO-IONOTENS using the switch on the right side of the device;

## STORED PROGRAMS

Read the following instructions if you want to use stored programs on MIO-IONOTENS.

### 1. PROGRAM SELECTION

Press MODE to scroll up/down programs list from P01 to P16, and stop on choosed program (see "Programs list" paragraph for details);

### 2. INTENSITY REGULATION

Use the two knobs on the sides of the unit to adjust intensity level for channel 1 and channel 2, up to 100mA max. For each channel (40mA max for ionophoresis programs). The intensity value can be adjusted in 1mA each step.

For patient safety, if the cables are not correctly connected MIO-IONOTENS resets intensity value

when it reaches 10mA and the display will turn to red light.

Display shows countdown timer till the end of the program. An output signal advice of the end of the program.

## USER PROGRAMS

You can follow the following steps to set a new program with your own preferred values or the ones recommended by your doctor/physiotherapist.

### 1. CHOOSE PROGRAM

Press MODE to scroll up/down programs list and choose one of the four user program P1, P2, P3, P4, then pres SET to enter menu and set the parameters for the choosed program;

### 2. SET PARAMETERS

- Adjust waveform frequency (0,5÷150 Hz), increasing or decreasing value using ▲(increase) and ▼(decrease);
- Press SET to confirm;
- Adjust impulse amplitude (30÷300µs), increasing or decreasing value using ▲(increase) and ▼(decrease);
- Press SET to confirm;
- Adjust therapy time (1÷60 minutes or continuous), increasing or decreasing value using ▲(increase) and ▼(decrease);
- Press SET to confirm.

### 3. INTENSITY REGULATION

Use the two knobs on the sides of the unit to adjust intensity level for channel 1 and channel 2, up to 100mA max. For each channel (40mA max for ionophoresis programs). The intensity value can be adjusted in 1mA each step.

For patient safety, if the cables are not correctly connected MIO-IONOTENS resets intensity value when it reaches 10mA and the display will turn to red light..

Display shows countdown timer till the end of the program. An output signal advice of the end of the program.

Press MODE during treatment to exit therapy: MIO-IONOTENS will show the programs list screen.

Press SET during treatment to show the program in use.

WARNING. Press SET and ▲ buttons at the same time to lock the setting parameters: a key icon will appear on the display. Press them again to unlock.


WARNING. MIO-IONOTENS will enter into “sleep mode” if not press any key for 2 minutes to preserve batteries.

WARNING. To preserve batteries we suggest to use ionophoresis programs (P01, P02, P03, P15 and P16) only on channel 1.

## TENS and ionophoresis

For TENS programs, the intensity should be adjusted to a level between the thresholds of perception and pain. The maximum limit is reached when the muscles surrounding the area treated begin to contract. It is best to stay below that limit.

The electrodes should be positioned to form a square surrounding the painful area using Channel 1 and Channel 2 as shown in illustration 1.

For ionophoresis programs set up an intensity value so to have “pins and needles” on treatment area.

The used drug can have negative polarity, positive polarity or double polarity. The current induce the drug to run from one pole to the other, crossing the painful area and releasing the specific active ingredient.


WARNING: before starting the therapy, wet the sponge electrodes and wring them so to avoid dripping, then put the drug on the electrode as follow:

- Active polarity drug: dissolve this drug on the electrode connected to positive pole (red connector)
- Negative polarity drug: dissolve this drug on the electrode connected to negative pole (black connector)
- Double polarity drug: can be dissolved on positive pole or negative pole without distinction.

Position the electrode with the drug on painful area, and the other electrode on the other side.

At the end of the program, the skin could lightly turn bright red; the reddening usually vanishing few minutes after the end of program.

WARNING. Do not use ionophoresis program near metallic prosthesis.

### LIST of drugs used in ionophoresis

Drug	Polarity	Pharmaceutical action	Indications
Calcium chloride (Sol, 1% 2%)	Positive	Sedative and recalcifying	Osteoporosis, Spasmophilia, algodystrophic syndrome DO not use in cases of arteriosclerosis
Magnesium chloride (Sol. 10%)	Positive	Analgesic, sedative, Fibrolitic	Substitute for calcium chloride cases with arteriosclerosis
Potassium iodide	Negative	Sclerolytic, emollient	Scars, adhesions, Dupuytren's disease, cheloidis
Lysine acetylsalicylate	Negative	Antiphlogistic and analgesic	Artrosis
Flectadol, Aspegic	Negative	Analgesic	Artrosis extra/intra- articular rheumatism
Local anaesthetics (Novocaine, lidocaine)	Negative		Local anaesthesia, trigeminal neuralgia
Benzydamine	Positive	Antiphlogistic and analgesic	Rheumatoid arthritis

Diclofenac sodium	Positive/Negative	Antiphlogistic and analgesic	Bruising
Orudis, voltaren, Feldene, Lometacen, Arfen, Tilcotil, Axera, Naprosyn	Negative	Anti-inflammatory	Degenerative and extra-articular rheumatism, gout
Piroxicam	Positive	Antiphlogistic and analgesic	Sprains
Sodium salicylate (1%-3%)	Negative	Analgesic	Articular rheumatism, acute decongestant, myalgia
Ketoprofene Lysine salt	Positive/negative	Anti-inflammatory	Arthrosis, arthritis
Thiomucase	Negative	Antiedemic	Post-trauma and post surgical edema due to venous insufficiency

If the drug used is not included in the above list, determine the polarity from the package or consult the prescribing doctor or dispensing pharmacist.

### Stored programs

PRG	Hz	Impulse width	Therapy time	Program details
P1	0.5-5Hz adjustable	30-300µs	1-60 min o continuous	Burst impulse
P05	2Hz	180µs	30 min	Burst impulse
<p>This program produces a training effect using the frequencies of conventional TENS. Useful for pain therapy. The action is similar to the one of endorphinic TENS, with a treatment time of 15 minutes for a phase.  <b>Position of electrodes:</b> form a square above the painful area as shown in illustration 1.</p>				

PRG	Hz	Impulse width	Therapy time	Program details
P2	1-150Hz adjustable	30-300µs	1-60 min o continuous	Conventional TENS
P04	80Hz	180µs	30 min	Conventional TENS
P11	80Hz	60µs	30 min	Conventional TENS
P14	100Hz	75µs	30 min	Conventional TENS
<p>Its purpose is to induce the organism into blocking pain at the spine, in accordance with the "Gate Control Theory" by Melzack and Wall. Conventional tens activates large-diameter nerve fibres, blocking the path of small-diameter nerve fibres at the spine. So action is mainly taken against the symptom: to simplify it further, the wire transmitting pain information is obstructed. Treatment duration should be no less than 30/40 minutes. Conventional tens is a current that can be used to treat general everyday pain. The average number required to benefit from the treatment is 10/12 per day (there are no contraindications for up to double this amount).  <b>Position of electrodes:</b> form a square above the painful area as shown in illustration 1.</p>				

PRG	Hz	Impulse width	Therapy time	Program details
P3	1-150Hz adjustable	30-300µs	1-60 min o continuous	Modulated pulse width on 10 seconds period, from settled width value to 60% less
P06	80Hz	50-180µs	30 min	Modulated pulse width
<p>Modulated pulse width stimulation is a type of high-frequency stimulation where the pulse width varies continuously. This can cause an undulating sensation, which may be more pleasant than a constant pulse duration. Program intended for pain alleviation and a massage effect on muscle like the trapezius.</p> <p><b>Position of electrodes:</b> form a square above the painful area as shown in illustration 1.</p>				

PRG	Hz	Impulse width	Therapy time	Program details
P4	1-150Hz adjustable	30-300µs	1-60 min o continuous	10 seconds modulated frequency, from settled frequency value to 40% less
P13	2-80Hz	200-100µs	30 min	Frequency modulation and amplitude width modulation
<p>Modulated frequency stimulation is a type of high-frequency stimulation where the frequency varies continuously. This can cause an undulating sensation, which may be more pleasant than a constant pulse duration. Program intended for pain alleviation and a massage effect on muscle like the trapezius.</p> <p><b>Position of electrodes:</b> form a square above the painful area as shown in illustration 1.</p>				

PRG	Hz	Impulse width	Therapy time	Program details
P01	800Hz	100µs	20 min	Standard ionophoresis 1, monophasic impulse
P02	1000Hz	100µs	20 min	Standard ionophoresis 1, monophasic impulse
P03	1200Hz	100µs	20 min	Standard ionophoresis 1, monophasic impulse
P15	1100Hz	100µs	20 min	High impact ionophoresis 1
P16	1300Hz	100µs	20 min	High impact ionophoresis 2
<p>For ionophoresis program the intensity must be strong enough to produce a relevant perception, near pain, till the muscles surrounding the area treated begin to contract.</p> <p><b>Position of electrodes:</b> position the electrode with the drug on painful area, and the other electrode on the other side.</p>				

PRG	Hz	Impulse width	Therapy time	Program details
P07	2Hzx3sec + 15Hzx3sec	180µs	30 min	Endorphinic Tens + slow twitch stimulation
P08	2Hzx3sec + 80Hzx3sec	180µs	30 min	Endorphinic Tens + fast twitch stimulation
P09	100Hzx3sec +	150µs+200µs	30 min	Waveform, amplitude width and

	2Hzx3sec			frequency modulation
P12	2Hz	250µs	30 min	Endorphinic TENS
<p>This type of stimulation produces two types of effects according to how the electrodes are positioned: positioning the electrodes in the dorsal region, promotes the endogenous production of morphine-like substances capable of raising the pain perception threshold; positioning the electrodes to form a square above the painful area as shown in illustration 1 produces a vascularizing effect. Vascularization increases arterial flow and consequently aids the removal of algogenic substances and helps to restore normal physiological conditions. Do not position the electrodes close to inflamed areas. Intensity adjusted for good solicitation of the part stimulated, the sensation must be similar to that of a massage. <b>Position of electrodes:</b> form a square above the painful area as shown in illustration 1.</p>				

PRG	Hz	Impulse width	Therapy time	Program details
P10	10Hz	180µs	30 min	Nausea treatment
<p>Low frequency TENS similar to endorphinic TENS, particularly indicated for nausea treatment. <b>Position of electrodes:</b> place the black cable electrode inside of the wrist and the red one outside of the wrist.</p>				

### 25 Cervicale

(per programma medicale)

- Ⓢ Area: cervicale (for medical programme)
- Ⓣ Zone: cervicale (pour programme médical)
- Ⓛ Bereich: Halswirbel (medizinisches Programm)
- Ⓧ Zona: cervicale (para programa médico)
- Ⓝ Zone: Hals (geneskundig programma)

35 x 45 mm  
45 x 50 mm

35 x 45 mm  
45 x 50 mm

### Spalla (periartrose) 26

(per programma medicale)

- Ⓢ Area: shoulder (peri-arthritis) (for medical programme)
- Ⓣ Zone: épaule (périarthrite) (pour programme médical)
- Ⓛ Bereich: Schulter (Periarthritis) (medizinisches Programm)
- Ⓧ Zona: hombro (periartrosis) (para programa médico)
- Ⓝ Zone: schouder (periartrosis) (geneskundig programma)

35 x 45 mm  
45 x 50 mm


35 x 45 mm  
45 x 50 mm

### Lombare (lombalgia) 27

(per programma medicale)

- Ⓢ Area: lumbar (lombalgia) (for medical programme)
- Ⓣ Zone: lombaire (lombalgie) (pour programme médical)
- Ⓛ Bereich: Lenden (Lendenschmerzen) (medizinisches Programm)
- Ⓧ Zona: lumbar (lumbago) (para programa médico)
- Ⓝ Zone: lenden (lage rugpijn) (geneskundig programma)

35 x 45 mm  
45 x 50 mm

35 x 45 mm  
45 x 50 mm

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(per programma medicale)

- Ⓢ Area: the length of the leg (sciatica) (for medical programme)
- Ⓣ Zone: la longueur de la jambe (sciatalgie), (pour programme médical)
- Ⓛ Bereich: Bein in der Länge (Ischias) (medizinisches Programm)
- Ⓧ Zona: la longitud de la pierna(cíatica) (para programa médico)
- Ⓝ Zone: het hele been (ischialgie) (geneskundig programma)

35 x 45 mm  
45 x 50 mm

35 x 45 mm  
45 x 50 mm

**29 Gomito (epicondilitis)**  
(per programma medicale)

Ⓞ Area: elbow (epicondylitis)  
(for medical programme)

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Ⓞ Area: hip (hip pain)  
(for medical programme)

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Ⓞ Area: wrist (carpal tunnel)  
(for medical programme)

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**Maintenance**

**Replacing the batteries**

In case of low batteries, display shows a flashing icon  on the right lower side: take care that could be not possible to start the therapy or complete it.

For batteries replacement please follow the steps below:

- Switch off MIO-IONOTENS using ON/OFF switch;
- Open batteries compartment on the back side of MIO-INOTENS and remove the batteries;
- Insert 4 new 1,5V AAA Size Alkaline batteries;
- Close batteries compartment.

WARNING. MIO-IONOTENS can work with Ni-MH 1,2V AAA Size rechargeable batteries: follow manufacturer instructions for the correct use.

WARNING. for disposal exhausted battery follow the instructions on chapter “Disposal”. Do not open or burn battery. Do not short circuit poles. Keep battery away from sparks or naked flames. In the event of internal electrolyte coming into contact with skin or garments, wash immediately with water. In the event of electrolyte coming into contact with eyes, rinse thoroughly and seek medical attention.

## Cleaning

Clean the equipment from the dust using a dry soft cloth.

Resistant stains can be removed using a sponge soaked in solution of water and alcohol.

## Carriage and storage

### Carriage precautions

MIO-IONOTENS is a portable device, so it does not need any particular carriage precautions.

However we recommend to put away MIO-IONOTENS and its accessories in their own bag after every treatment.

### Storage precautions

MIO-IONOTENS is protected till following environmental conditions:

Outside of the packaging

Temperature	from +5 to + 40 °C
Humidity	from 30 to 80%
Pressure	from 500 to 1060 hPa

Inside of the packaging

Temperature	from -10 to +55 °C
Humidity	from 10 to 90%
Pressure	from 500 to 1060 hPa

## Disposal

The equipment is subjected to WEEE regulations (see the symbol  on the label) concerning separate waste collection: when disposing this product, please use the designed areas for disposing electronic waste or contact the manufacturer.

## Troubleshooting

If it is used in accordance with the instructions of the user manual, MIO-IONOTENS does not need a particular regular maintenance.

If you find any malfunctioning using MIO-IONOTENS, please follow these instructions:

- **MIO-IONOTENS doesn't switch on and/or the display doesn't light on.** Check batteries status and eventually replace them (see paragraph "replacing the batteries"). If the problem persists contact the manufacturer.
- **MIO-IONOTENS doesn't transmit electric impulses.** Check that the cable jacks have been inserted in the electrodes and that the plastic protection has been removed from the electrode. Check that the cables have been connected correctly (connector well inserted in the device). Check cables and electrodes status (no signs of damage) . If the problem persists contact the manufacturer.
- **MIO-IONOTENS switches off during operation.** Replace the batteries and start the treatment again. If the problem persist contact the manufacturer.
- **MIO-IONOTENS doesn't permit to increase output intensity or does not keep the settled value and resets.** Replace the batteries and start the treatment again. If the problem persists contact the manufacturer.

## Repairing

Every intervention on device must be performed by manufacturer or authorized personnel.

## Spare parts

Contact the manufacturer or local distributor for original spare parts.

To preserve product warranty, functionality and product safety we recommend to use only original spare parts.

**Warranty (we recommend that you fill in and post the light blue warranty card).**

MIO-IONOTENS is covered by a 24-month warranty with effect from the date of purchase. Connection cables are covered by a 6-month warranty; parts subject to wear (like electrodes and batteries) are excluded from the warranty, unless evident manufacturing defects are found.

The warranty will lapse in the event that the equipment is tampered with or work has been carried out on the same by personnel not authorised by the manufacturer or authorised dealer.

The warranty conditions are described in the “Warranty Regulations”.

### **Post the light blue card and keep the green card**

N.B. responsibility of the purchaser: the purchaser is responsible for validating the warranty by filling in the warranty card (light blue warranty validation card), stamped by the retailer, and sending it to the manufacturer.

In the event of future repairs under warranty, the equipment must be packaged to prevent damage during transport and sent to the manufacturer together with all of the accessories, ensuring that all parts of the service card (green service card) have been filled in. The purchaser only has the right to repair under warranty when the equipment is returned to the manufacturer complete with the green service card signed by the retailer, together with the receipt or invoice proving the correct origin of the product.

### **Warranty regulations**

1. We recommend that the purchaser fills in the light blue warranty validation card, stamped by the retailer, and posts it to the manufacturer immediately after purchase.
2. In the event of repairs under warranty, the purchaser must fill in all parts of the green service card and send it to the manufacturer together with the equipment. The receipt for tax purposes or the purchase invoice must be sent together with the goods.
3. The electronic parts are covered by a 24-month warranty. The warranty is given through the point of sale or directly from the manufacturer.
4. The warranty covers exclusively product damage causing operational defects.
5. The warranty covers exclusively the repair or replacement free of charge, including labour, of components found to be defective in terms of manufacture or material.
6. The warranty does not apply to damage caused by neglect or use not complying with the instructions provided, damage caused by work carried out by unauthorised personnel or damage caused by accidental causes or the buyer’s negligence, with particular reference to external parts.
7. The warranty does also not apply to damage to the equipment caused by incompatible power supplies (the equipment operates at 4.8V with an internal battery pack).
8. Parts subject to wear after use are excluded from the warranty.
9. The warranty does not include transport costs to be paid by the purchaser in relation to the method and speed of transport.
10. The warranty expires after 24 months. After such time repair work will be carried out at the rates currently in force for the parts replaced and the labour and transport costs.
11. Any controversy will fall within the exclusive jurisdiction of the Venice courts.

**I.A.C.E.R S.r.l. SU**

*Sede operativa:*

30030 Martellago (VE) - Via. S. Pertini 24/A  
Tel +39 041 5401356 - Fax +39 041 5402684

-----  
*Sede legale:*

S. Marco 2757 - 30124 Venezia  
Cod. Fisc./P.IVA IT 00185480274  
R.E.A. VE N. 120250 - M. VE001767  
Cap.Soc. € 110.000,00 i.v.  
w w w . i a c e r . v e . i t - i a c e r @ i a c e r . i t