GIMA

TERMOMETRO AURICOLARE TH889F
TH889F EAR-TYPE THERMOMETER
OHRENTHERMOMETER TH899F
THERMOMETRE AURICOLAIRE TH889F
TERMÓMETRO AURICULAR TH889F
TERMÔMETRO AURICULAR TH889F
ΘΕΡΜΟΜΕΤΡΟ ΑΥΤΙΟΥ TH889F

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
GEBRAUCHS- UND WARTUNGSHANDBUCH
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ

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
PREMISE

Thank you for purchasing this GIMA Ear-type thermometer. This thermometer has been planned with advanced infrared technologies and room environment compensation, for accurate and prompt measurements.

Do not use this device for Live & Die decision or Safety related applications.

FUNCTIONING

Before taking any measurement, make sure that the ear canal is clean, dry and not obstructed.

1) Pull off the probe cap **(1)**. Do not force in removing the cap.

2) It is necessary to apply a disposable plastic cap **(2)** on the probe **(3)** by means of the apposite connection ring **(4)**:

- Place the plastic cap **(2)** on the connection ring **(4)** in conformity with the direction indicated in the picture;

- Push the ring on the probe until a click sounds;

- Once the cap has been fixed, delicately remove the connection ring **(4)**.

Always make use of a new disposable plastic cap and make sure that it has not been damaged.

3) Press ON/MEM **(6)** pushbutton. When the ear icon stops flashing and the equipment beeps twice, the thermometer is ready for use.

4) Gently pull the ear back to straighten the ear canal and snugly position the probe **(3)** in the ear canal, aiming towards the membrane of the eardrum to obtain an accurate reading.

5) Measuring the ear temperature:

Hold the SCAN **(5)** pushbutton pressed with your forefinger until you hear a beep. At this point the LCD display **(7)** shows the detected temperature.

5) Device will automatically power off if left idle for more than 60 seconds to extend battery life.

SWITCHING BETWEEN FAHRENHEIT AND CELSIUS

To change the unit of measurement of temperature from °C to °F, hold down the SCAN button **(5)**, when the unit is off, then press and hold the ON/MEM button **(6)** pressed until the display shows the symbol °F or °C. Then release both buttons. Use the same process to change the unit of measurement again.

MEMORY FUNCTION

If the reading thermometer is within the normal temperature range ($34^{\circ}\text{C} \div 42,2^{\circ}\text{C}$ or $93,2^{\circ}\text{F} \div 108^{\circ}\text{F}$), the data is saved into memory each time the thermometer is turned off.

Press the ON/MEM button (6) when the thermometer is turned on to see the temperature stored.









WARNING

- Follow the instructions of the manual to the letter for the best measurement.
- It is advisable to stay in the same environment for at least 5 minutes before taking the temperature.
- Do not have a bath or physical exercise before a measurement.
- While measuring avoid exposing to direct sunlight, rain or wind; do not take measurements near source of heat or very humid places.
- Keep far from children's reach both the connection ring (4) and the disposable plastic caps (2).

BATTERY REPLACEMENT

This device is supplied with one lithium cell (CR2032 x 1). Remember to power off the unit before replacing the battery. A malfunction may occur if the power is on when the battery is replaced. If a malfunction occurs, restart the device. Keep the battery away from children.

PROBLEM SOLVING

Messaggio d'errore	Problema	Soluzione
	Device stabilization in process.	Wait until  stops flashing.
	Battery is low and no more measurements are possible.	Replace the battery.
	Measurement before device stabilization.	Wait until  stops flashing.
	The device showing a rapid ambient temperature change.	Allow the thermometer to rest in a room at least 30 minutes at room temp. ($10^{\circ}\text{C} \div 40^{\circ}\text{C}$, $50^{\circ}\text{F} \div 104^{\circ}\text{F}$).
	The ambient temperature is not within the range between 10°C and 40°C ($50^{\circ}\text{F} \div 104^{\circ}\text{F}$).	Allow the thermometer to rest in a room at least 30 minutes at room temp. ($10^{\circ}\text{C} \div 40^{\circ}\text{C}$, $50^{\circ}\text{F} \div 104^{\circ}\text{F}$).
	Error 5 ÷ 9, the system is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.
	Temperature taken is higher than $42,2^{\circ}\text{C}$ (108°F).	Check the integrity of the probe cover and take a new temperature measurement.
	Temperature taken is lower than 34°C ($93,2^{\circ}\text{F}$).	Make sure the probe cover and the lens are clean and take a new temperature measurement.
	Device can not be powered on to the ready stage.	Change with a new battery.

CLEANING AND STORAGE

Replace the disposable plastic cap after each use to ensure an accurate reading and avoid cross infections.

In case multiple consecutive measurements are taken on the same patient, clean the probe plastic cap (2) with a soft cloth moistened with medical alcohol.

The probe lens (8) is the most delicate part of the thermometer. Use with care during cleaning to avoid damage. If the device is accidentally used without a plastic cap (2), clean the probe lens as:

- a) Gently wipe lens surface with a cotton swab or soft cloth moistened with water or medical alcohol;
- b) Allow the lens to fully dry at least 30 minutes before attaching a new probe cap.

Always store at a temperature between 10°C to 40°C (50°F ÷ 104°F), in a dry place and far from liquids and direct sunlight. In case extraordinary maintenance is required, apply to the closest retailer.

GUARANTEE

Thank you for having acquired our product. This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, doppler probes, SpO₂ sensors, lamps, batteries, electrodes, handpieces, washers and packings, resistances, etc.

No compensation will be given for loss of use of the product. Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.

The guarantee will be invalidated if the registration number has been removed, cancelled or altered.

Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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INNOVATION
RADIANT

RADIANT INNOVATION INC.
<http://www.radiantek.com.tw>

DECLARATION OF CONFORMITY(Annex VII and V, Class Iia products; with product list)

Radiant Innovation Ear Thermometer and Probe Cover

Manufacturer: RADIANT INNOVATION INC.

*Address: 1F, No.3, Industrial E. 9th Rd., Science-Based Industrial Park,
HsinChu, Taiwan 300, R.O.C.*

DECLARATION OF CONFORMITY
medical devices

We hereby declare that the distributed CE marked products, specified in the annexed product list, are covered by the "CE Marking of Conformity Certificate", reference number:HD 60005278 0001, expiry on Jun 1, 2006 and delivered by TUV Rheinland Taiwan, Notified Body Identification Number 0197, and conform to the required technical documentation, in accordance with Annex VII of the "EC-Directive", the Council Directive 93/42/EEC of 14 June 1993, concerning medical devices.

In addition, we ensure and declare that the distributed CE marked products, as mentioned and falling within Class IIa, meet the provisions of the EC-Directive which apply to them.

This declaration is based on the application of the Quality System approved for the manufacture and final inspection of the products concerned, in accordance with Annex V of the EC-Directive. The conformity of the production quality assurance set out in Annex V, is described in the said CE Marking of Conformity Certificate, issued and delivered by TUV Rheinland.

This declaration is supported by the Quality System certification based on the harmonized standards ISO 9001:2000 / ISO 13485:2000, Quality System Certificate with reference number: SY 60005279 0001, expiry on Jun 1, 2006 and delivered by TUV Rheinland.

This Declaration of Conformity covers *Radiant Innovation* as specified in the product list belonging to this declaration, and is valid for all products concerned bearing the CE marking and manufactured at the following site(s):

RADIANT INNOVATION INC.

**Factory: 1F, No.3, Industrial E. 9th Rd., Science-Based Industrial Park,
HsinChu, Taiwan 300, R.O.C.**

**Factory: No.468, XiaoLin Road, ChengBei, KunShan City, JianSu, China.
TEL:886+3+564 4185 FAX:886+3+564 4170**

Issued on Jul.22, 2003
Hsin-Chu, Taiwan R.O.C.

Jolin Yen / President
RADIANT INNOVATION INC.

1F, No.3, Industrial E. 9th Rd., Science-Based Industrial Park, HsinChu, Taiwan 300, R.O.C.

● TEL: 886+3+5644185 ● FAX: 886+3+5644170