

**Nome Prodotto**

Imbracature Deluxe 300 con Testiera in PVC.

Codice	Misure	Portata (Kg)
01.01.2036.300	S	300
01.01.2037.300	M	300
01.01.2038.300	L	300
01.01.2039.300	XL	300

**Destinazione d'uso**





Le Imbracature vengono associate ad un sollevatore della Linea SCANMEDICALEUROPE per la movimentazione dei pazienti.

**Descrizione**

L' IMBRACATURA DELUXE 200 CON TESTIERA, grazie al suo design, consente il trasferimento del paziente da seduto a sdraiato e viceversa.

L' appoggiatesta fornisce un comodo supporto per la nuca prevenendo movimenti inappropriati della cervicale.

Per garantire la massima sicurezza, le imbracature dovrebbero essere sostituite ogni 3 anni.

Materiali	
Imbracatura	PVC
Manutenzione e Pulizia	
	Lavaggio Max. 40°C
	Non Stirare
	Non Candeggiare
	Non Asciugare a Tamburo
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.a.S.
Classe di Dispositivo 93/42	Classe I
Codice ISO	12.36.21.006
Codice CND	Y123603
Codice RDM	1301118/R

**Product Name**

PVC Deluxe 300 with Head Support.



Code	Size	Capacity (Kg)
01.01.2036.300	S	300
01.01.2037.300	M	300
01.01.2038.300	L	300
01.01.2039.300	XL	300

**Use**





The slings are associated with a lift of the SCANMEDICALEUROPE Line for patient handling.

**Description**

Its design allows transfers of patients who are in a sitting or lying position.

The head support prevents inappropriate cervical movements.

To ensure maximum safety, the sling should be replaced every 3 years.

Materials	
Sling	PVC
Maintenance and Cleaning	
	Water Temperature not Above 40°C or 105°F
	Do not Iron
	Do not Bleach
	Do not Tumble Dry
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.a.S.
Directive 93/42	Class I
ISO Code	12.36.21.006
CND Code	Y123603
RDM Code	1301118/R